

Vernacare™

solo

The New Sleek and Compact Single-Item Macerator

A WORLD LEADER IN HOSPITAL HUMAN WASTE MANAGEMENT



The Solo macerator is ideal for private rooms providing discreet and efficient elimination of biodegradable utensils and human waste.

Solo's simple, streamline and aesthetically appealing appearance makes it perfect for integration into the décor of private rooms. The easy-to-use Solo reliably macerates pulp units quickly and quietly with the press of a button. The in-room Solo helps to control unpleasant odours, by immediately eliminating waste which also helps to reduce cross contamination and enhance infection control.

Solo offers many benefits that make it an ideal in-room macerator:

- Fits in a patient's washroom eliminating the need to cross the infection control barrier
- Suitable for low volume usage in private rooms, isolation rooms, critical care rooms, outpatient clinics and offsite dialysis centres
- Flexible and ideal for long term care facilities
- Easy to use, install and maintain
- Helps reduce cross contamination and enhances infection control
- Energy efficient by using cold water and quick cycle times
- Quiet motor is ideal for night time use

Solo Product Information Details

Item Number:	V40-400
Electrical Requirements:	110 volts
Cycle Time:	Less than 50 seconds
Water Requirements:	Cold water only
Capacity:	One Vernacare pulp utensil

For more product specification information, visit www.vernacare.com



Vernacare products are biodegradable and made from 100% recycled post-consumer newsprint.


Vernacare™

solo

The New Sleek and Compact Single-Item Macerator



Solo Product Dimensions




Installation Notes

1. Drain outlet terminates at a fitting suitable for a 38 mm (1 1/2") dia. pipe stub. P trap supplied and fitted internally.
2. Water supply in via 700 mm (27.56") long 10 mm (0.394") bore flexible water supply hose with 3/4" GHT female thread supplied.
3. Electrical Supply: Connect to a 15 Amp Circuit with Motor Switch.

Warranty: 1 Year, Parts Only.

For more product specification information, visit www.vernacare.com



CONTACT YOUR VERNACARE REPRESENTATIVE FOR MORE INFORMATION.

Vernacare Head Office
150 Norfinch Drive, Unit 4
Toronto, Ontario
Canada M3N 1X6

Toll Free: 1-800-268-2422
Toll Free Fax: 1-877-311-2811
Telephone: 416-661-5552
Fax: 416-661-5559

Email: marketing@vernacare.com
Website: www.vernacare.com

ENVIRONMENTAL PRINTING – IT'S MORE THAN RECYCLED PAPER.™

